



A.R.I.D.

STAFF OPERATIONS MANUAL

Web Dashboard User Guide

A practical guide for LGU personnel and health workers who monitor, assess, and review synchronized dengue-vector breeding-site reports.

Platform	Web dashboard
Audience	Authorized staff users
Edition	August 2026

About this guide

This guide documents the web dashboard currently available in A.R.I.D. It explains the staff workflow from signing in through reviewing reports and analyzing separate site photos.

The web dashboard is for monitoring. Field capture, GPS tagging, classification, and report submission happen in the A.R.I.D. mobile app. New reports appear on the dashboard after they synchronize.

Contents

01 Access and sign in

02 Monitor reports on the map

03 Inspect and review a report

04 Analyze a separate photo

05 Status messages and troubleshooting

06 Daily-use checklist

Internet connection

Required for sign-in, live reports, report photos, and saving review-status changes.

Authorized account

Your signed-in account must be included in the dashboard staff list.

Access and sign in

1 Open the web dashboard

Use a current desktop or mobile browser. Confirm that the device has an internet connection.

2 Enter your staff credentials

Type your work email and password, then select **Sign in**. If your organization enabled Google authentication, you may instead select **Continue with Google**.

3 Wait for access verification

The dashboard checks both your sign-in and your staff authorization before showing monitoring data.

If **"Access not enabled yet"** appears: your sign-in succeeded, but your account is not yet on the staff list. Copy the account UID shown on screen and send it to the dashboard administrator. After access is added, select **Check again**; you do not need to sign out.

Useful account actions

ACTION	WHEN TO USE IT
Switch	You signed in with the wrong account and need to choose another one.
Copy UID	An administrator needs your unique account identifier to enable staff access.
Check again	An administrator has added your access and you want the page to verify it immediately.
Sign out	You have finished using the dashboard, especially on a shared device.

Monitor reports on the map

The monitor page combines scope filters, risk controls, map layers, a live map, an area snapshot, and a contributor leaderboard. New synchronized reports appear automatically.

Choose which reports are included

CONTROL	PURPOSE
Study area	Limits the view to reports inside the selected geographic boundary.
Date range	Shows all dates, the last 24 hours, the last 7 days, or the last 30 days.
Classification	Shows all reports, breeding-site results, or non-breeding results.
Risk	Shows all, high, moderate, or low-risk reports. The displayed number is the count for that choice.

Choose how reports are drawn

Markers

Individual report locations. Select a marker to inspect its report.

Heatmap

A density view of reports currently included by the filters.

Hotspots

Clustered report areas that may warrant closer attention.

Read the side panel

- **Area snapshot** shows the visible report total, risk distribution, recent activity, and weekly change.
- **Top contributors** ranks up to eight synchronized reporters by points.
- Counts and summaries update to match the selected area, date, classification, and risk filters.

Inspect and review a report

1 Select a report marker

Choose a marker on the map. The report detail opens at the top of the right-side panel.

2 Assess the available evidence

Review the field photo, risk level, breeding/non-breeding classification, model confidence, GPS coordinates, GPS accuracy, reporter, points awarded, and synchronization state.

3 Set the review status

Select **Unreviewed**, **Reviewed**, or **Actioned** according to the work completed.

4 Close the report or continue filtering

Use the close button to return to the summary. If a filter hides the selected report, the dashboard warns you but keeps the detail open until you close it.

Unreviewed

The report is awaiting staff assessment.

Reviewed

A staff member has checked the report evidence.

Actioned

A response or follow-up has been carried out.

Data protection: review status is the only report field changed by the dashboard. The classification and risk level remain as recorded by the mobile app.

Analyze a separate photo

Use **Analyze image** when you want to inspect a photo that is not part of a synchronized report.

1 Open image analysis

Select **Analyze image** in the monitor header.

2 Choose a supported file

Drop a photo into the upload area or select **Choose image**. JPG, PNG, and WebP images up to 15 MB are supported.

3 Wait for local analysis

The first analysis also loads a 27 MB model. The selected photo stays in the browser and is not uploaded.

4 Interpret the result

Detected objects are outlined on the photo and listed with confidence percentages. The model recognizes bottles, coconut exocarps, drain inlets, tires, and vases.

5 Verify the site

Use the result to support prioritization, then confirm the actual site conditions through field inspection.

A detection is not confirmation of a breeding site. It indicates a potential breeding container only. Confirm stagnant water, larvae, or mosquitoes in the field.

No detection is not proof of no risk. Objects below the 25% confidence threshold are not displayed, and the model may not recognize every relevant condition.

Status messages

Recently updated

The dashboard is connected and live report data is current.

Offline

The map shows the last data received. Reconnect to receive new reports or save review changes.

Live updates stopped

The existing view remains visible, but new data is paused. Check the connection and follow the message displayed.

Photo unavailable

The report has no photo URL, or its Cloudinary image could not be loaded. Use the remaining report details and retry later.

No report selected

Select a marker to inspect its photo, confidence, GPS information, and review status.

No target found

Nothing exceeded the detector's confidence threshold. This does not establish that the location is risk-free.

Appearance: use the sun/moon button to switch between light and dark display modes. The choice does not change report data.

Quick staff checklist

Use this short checklist at the beginning and end of each monitoring session.

- ✓ Confirm that the dashboard is online and recently updated.

- ✓ Select the correct study area and date range.

- ✓ Prioritize high-risk reports and visible hotspots.

- ✓ Open report markers and verify the photo, confidence score, and GPS quality.

- ✓ Set the appropriate review status after assessment or response.

- ✓ Treat browser detections as screening results and verify conditions in the field.

- ✓ Sign out before leaving a shared workstation.

Remember: the dashboard supports operational judgment; it does not replace field verification or local dengue-control protocols.

A.R.I.D.

Autonomous Risk Identification and Distribution Mapping

Web dashboard documentation for authorized LGU and health-worker staff.